

CHILI PACS

IHE Integration Statement

Version 2019





Manufacturer
NEXUS / CHILI GmbH

Friedrich-Ebert-Str. 2
D-69221
Dossenheim/Heidelberg

Phone +49 6221 18079-10
Fax +49 6221 18079-11
E-mail info@nexus-chili.com

Internet: www.nexus-chili.com

Support

- + +49 6221 18079-20
- + support@nexus-chili.com

Revised: April 2019
CHILI PACS Version 4.64
IHE Integration Statement Version 2019



© 2003-2020 NEXUS / CHILI GmbH

The authors reserve the right to make changes at any time without prior notice. All rights reserved, in particular the right to reproduction, distribution, and translation. No part of the software or the user manual may be reproduced or processed, copied, or distributed by electronic means without the prior permission in writing from the authors.

IHE Integration Statement		Date: 2019-04-17
Vendor	Product Name	Version
NEXUS / CHILI GmbH	CHILI PACS	4.64
This product implements all transactions required in the IHE Technical Framework to support the IHE Integrations Profiles, Actors, and Options listed below:		

Integration Profiles Implemented	Actors Implemented	Options Implemented
Access to Radiology Information (ARI)	Report Repository	None
	External Report Repository Access	None
Audit Trail and Node Authentication (ATNA)	Secure Node	None
	Secure Node	Radiology Audit Trail
	Secure Node	FQDN Validation of Server Certificate
	Secure Node	BCP195 TLS Secure Transport Connection - All TLS Versions
	Secure Node	BCP195 TLS Secure Transport Connection - TLS 1.2 Floor
Basic Security (SEC)	Image Manager/Archive	None
Consistent Presentation of Images (CPI)	Print Composer	None
	Image Manager	None
	Image Display	None
Consistent Time (CT)	Time Client	None
	Time Server	None
Cross-Enterprise Document Media Interchange (XDM)	Portable Media Creator	CD-R, USB
	Portable Media Importer	CD-R, USB

Integration Profiles Implemented	Actors Implemented	Options Implemented
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Consumer	None
	Imaging Document Source	CDA Wrapped Text Report
	Imaging Document Source	PDF Report
	Imaging Document Source	Set of DICOM Instances
Evidence Documents (ED)	Image Display	None
	Image Manager/Archive	None
Healthcare Provider Directory (HPD)	Provider Information Directory	Provider Information Feed
Imaging Object Change Management (IOCM)	Image Display	None
	Image Manager/Archive	None
	Change Requester	None
Import Reconciliation Workflow (IRWF)	Importer	Scheduled Import - uses MWL, Unscheduled Import - no MWL
	Image Manager	None
Import Reconciliation Workflow (IRWF.b)	Importer	None
	Image Manager/Archive	None
Invoke Image Display (IID)	Image Display	None
Key Image Note (KIN)	Image Manager/Archive	None
	Evidence Creator	None
	Image Display	None
Mammography Image (MAMMO)	Image Manager	None
Patient Administration Management (PAM)	Patient Demographics Consumer	None
	Patient Demographics Consumer	Merge

Integration Profiles Implemented	Actors Implemented	Options Implemented
Patient Demographic Query (PDQ)	Patient Demographics Consumer	None
Patient Identifier Cross-referencing for MPI (PIX)	Patient Identifier Cross-reference Consumer	None
Patient Information Reconciliation (PIR)	Department System Scheduler/Order Filler	None
	Image Manager/Archive	None
Portable Data for Imaging (PDI)	Image Display	None
	Print Composer	None
	Portable Media Creator	None
	Portable Media Importer	None
Radiation Exposure Monitoring (REM)	Image Display	None
	Image Manager	None
Scheduled Workflow (SWF)	Department System Scheduler/Order Filler	None
	Image Display	None
	Image Manager/Archive	None
	Performed Procedure Step Manager	None
Scheduled Workflow (SWF.b)	Department System Scheduler/Order Filler	None
	Image Display	None
	Image Manager/Archive	None
	Performed Procedure Step Manager	None
Simple Image and Numeric Report (SINR)	Report Manager	None
Teaching File and Clinical Trial Export (TCE)	Export Manager	None
	Export Selector	None
Web-based Image Capture (WIC)	Image Capturer	JPEG Storage
	Image Manager	None

Integration Profiles Implemented	Actors Implemented	Options Implemented
Web-based Image Access (WIA)	Imaging Document Responder	None
	Imaging Document Source	None

Internet address for vendor's IHE information: http://www.chili-radiology.com/ihe		
Links to Standards Conformance Statements for the Implementation: http://www.chili-radiology.com/ihe/		
Links to general information on IHE		
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org